

Why use this foam over a liquid spray or wipes?

BevistoCryl is economical and safe. It is applied in controlled beads and maintains shape until wiped, rather than the overuse of liquid sprays which may evaporate or run off during interruptions, losing effectiveness and causing re-application. In some cases liquids have been known to overspray or run to the floor creating a slip hazard. Generally wipes are used very wastefully, being held incorrectly and disposed of before complete optimisation, whilst the supply gradually dries reducing the desired outcome.

Why is foam more economical than spray?

150ml of foaming solution is at a minimum, equivalent to 300ml of liquid spray. It is also more controlled and cost effective against solutions.

Is this a general all-purpose cleaner?

No. BevistoCryl a specialised foaming cleaner and disinfectant developed to address the need for a suitable solution to the clean and disinfect healthcare style mattresses, comfort chairs upholstery and vinyl type surfaces & for effective bacteria control with faster drying times, for powerful neutralising and re-odouriser, a ready to use solution for easy and safe use to assist in extending the life of sensitive healthcare surfaces.

BevistoCryl foam has a drawing action as well as the ability to break down soil and bio-burden on the surface, which then allows our blended biocides to get to work on the actual surface.

*Please note BevistoCryl can be used in spray bottles if you wish.