



Whiteness Perfect | #1

Carbamide peroxide at 10%, 16% and 22%.

**LEADERSHIP AND
HIGH ESTHETICS
LIKE NO OTHER**



**INTRODUCTORY
PRICE
25%
OFF!**



Desensitizing agents:

- Potassium Nitrate
- Fluoride



High level of water: keeps your teeth hydrated, reducing the risk of sensitivity.



Present in over 100 countries.



Safety: pH close to neutral.



Maximum effectiveness in 3 different concentrations: **10%, 16% and 22%.**



WEARING TIME

10% and 16%: 3 to 4 hours a day or overnight use.
22%: 1 hour a day.

“FGM’s Whiteness line is a reference in whitening gels. We can use Whiteness Perfect, for the take-home technique, in the concentration from 10% to 22%. I conducted several studies already published in articles and books with these materials and I use them daily in the office and in courses, normally in the concentrations of 10% and 16%.”

Prof. Dr. Fabiano Marson

Coordinator for the Specialization and Master program in dentistry, Faculdade Ingá, Maringá PR.




European Cosmetic certification European Union - Obelis.*

*concentrations 10% and 16%.

Regulation 1223 | 2009 | EC. Whitening products classified as medical device are not available in some countries of the European Community. Check product availability in your country.

Shade changes at different times of the evaluation for the whitening agent in the upper and lower arches.



Study proves that with the use of Whitening Perfect it is possible to achieve high rates of tooth whitening.

Source: Silva FMM, Nacano LG, Pizi ECG. Clinical Evaluation of Two Dental Whitening Systems. Rev Odontol Bras Central. 2012;21(57):473-9.

CLINICAL CASE

Authors: Prof. Dr. Amanda Vessoni Barbosa Kasuya, Prof. Dr. Patricia Saram Progiante, Prof. Dr. Juliana Delben, Prof. Dr. Cleveson de Oliveira e Silva and Prof. Dr. Fabiano Carlos Marson.



Initial



Final

AVAILABILITY

Kit 10% and 16%



Content: 5 syringes with 3g of gel each (15g of gel in the kit), 2 tray plates, tray case and applicator tips.

RRP \$75.00
NOW \$56.00

Indications:

10%
· Patients with mild to moderate pigmentation and higher propensity to dental hypersensitivity

Kit 22%



Content: 4 syringes with 3g of gel each (12g of gel in the kit), 2 tray plates, tray case and applicator tips.

RRP \$75.00
NOW \$56.00

16%
· Patients with moderate pigmentation, lower hypersensitivity and desire for faster tooth whitening.

Bulk



Content: 50 syringes with 3g of gel each, individually packed, in the 10%, 16% or 22% concentration.

RRP \$325.00
NOW \$245.00

22%
· Patients with intense pigmentation, lower hypersensitivity and a desire for faster tooth whitening.

Introductory offer valid for a limited time and not to be used in conjunction with any other offers.